

HUNTLEIGH



sonicaidTMFetalCare

Instructions for use

Document No. 775300-EN-7

Introduction	3
1.1 Intended use.....	3
1.2 Contra-indications	4
1.3 Cautions & warnings	4
2 Getting started	7
2.1 Screen layout	7
3 System operation	8
3.1 Home page.....	8
3.1.1 To return to the home page.....	8
3.1.2 Help	8
3.2 Bed & bed group functions	8
3.3 Adding a new patient:.....	10
3.4 Admit an existing patient	12
3.5 Multi-bed view	14
3.6 Reviewing stored patient data	15
3.7 Starting a CTG trace.....	16
3.8 Stopping a CTG trace	16
3.9 Single trace view	17
3.9.1 Trace details	18
3.9.2 Trace scrolling.....	19
3.9.3 Alerts	20
3.9.4 Adding a trace annotation.....	22
3.9.5 Printing a trace.....	23
3.9.6 Connecting SRX series Doppler - optional	24
3.9.7 CTG analysis (optional feature)	26
3.9.8 STAN.....	31
3.10 Switching views	33
3.11 Patient details view.....	34
3.11.1 Editing patient details.....	34
3.11.2 Retrieving stored CTGs.....	36
3.11.3 Retrieving stored analysis results	36
3.11.4 Adding patient notes.....	36
3.12 Chalkboard	38
3.13 Discharging a patient.....	39
3.14 Settings, Audit and Admin	39
4 Trouble-shooting	43
5 System maintenance	44
5.1 General maintenance	44
6 Service support	44
6.1 Licence renewal.....	44

Introduction

Sonicaid Fetalcare 3 is a full specification CardioTocoGraph (CTG) viewing & archiving system. Traces are automatically saved to the patient's record and are quick and easy to retrieve for case review.

The system provides the following key functions:

- Single bed live CTG view
- Multi-bed CTG view
- Patient management
- Automatic trace archiving
- Trace retrieval for review
- User configurable CTG alerts
- Optional CTG analysis
- Supports STAN fetal monitors
- Optional security and audit
- Optional Partogram
- Optional GDT
- Optional remote CTG
- chalkboard

1.1 Intended use

Sonicaid Fetalcare 3 is not a stand-alone medical device – it can only be used as part of a networked CTG system when connected to one or more fetal monitors.

It supports a range of makes/models of fetal monitor and there is no fixed limit to the number of beds/monitors that can be connected. Equally, there is no guarantee as to how many the system can support, as this will be determined by other factors such as server resource, network capacity and other local infrastructure limitations which are beyond our control.

Its primary intended uses include:

- Collection of real-time data from one or more fetal monitors.
- Displaying CTG traces.
- Archiving of CTG traces and associated data.
- Retrieval of archived traces and associated data for review.
- Alert functions relating to fetal parameters.
- Optional Dawes/Redman CTG analysis.
- STAN fetal monitors.
- Managing patients, entering notes, adding annotations and printing.

1.2 Contra-indications

Sonicaid Fetalcare 3 is not intended to provide a primary fetal monitoring function. Its use is secondary to the fetal monitor(s) and must not be relied on to alert users to potential problems, clinical issues or clinical management of the pregnancy. Clinicians retain full responsibility for all aspects of pregnancy management and for correct and effective use of the fetal monitor in the patient environment.

The alert functions are provided simply to alert users to FHR data which is outside user-set rate and time limits. They must not be relied on to alert users to problems like loss of contact or to clinical conditions such as tachycardia or bradycardia.

1.3 Cautions & warnings



3rd Party Software

Sonicaid Fetalcare 3 is designed to run on a dedicated PC/server. Under no circumstances must any 3rd party software be installed on the system without prior approval, in writing, from Huntleigh Healthcare. In the event of unauthorised software being installed, Huntleigh Healthcare cannot be held responsible for resulting data corruption, misrepresentation, loss or any other failing of the system.



While remote access clients may have 3rd party software installed, it is recommended that running multiple applications concurrently should be avoided, or kept to a minimum. If system resources are overloaded, data may be lost or corrupted. Huntleigh Healthcare cannot accept liability for any such loss or resultant problems or outcomes.



System back-up

As with any software system, crashes may occur at any time, and may result in loss or corruption of clinical data. Similarly, hardware failures may result in loss or corruption of data. While every effort is made to minimise this risk, it is strongly recommended that back-up protection measures are employed in line with industry standard practice. Huntleigh Healthcare cannot accept liability for any loss of data, corruption or other loss relating to data back-up failure or loss through any other causes.



Clinical management

Sonicaid Fetalcare 3 is not a diagnostic tool – it simply presents information. As with any computer / software system, bugs or faults may result in incorrect information being displayed. If any doubt as to fetal or maternal condition arises through using the Sonicaid Fetalcare 3 system, alternative measures must be undertaken immediately to ensure appropriate clinical management.



System Security

In the event of unauthorised access ('hacking') into the system or through any other malicious actions, data may be lost or corrupted. The system has no protection against unauthorised access. Users should take appropriate local measures to limit access to authorised users only.

**Data protection & patient confidentiality**

Due to the flexible, user configurable nature of Sonicaid Fetalcare 3, the system administrator is responsible for ensuring compliance with any local, national or other regulatory requirements relating to patient information, the storing, displaying and archiving of such data, and access to such data.

**Data integrity**

At all times, clinicians must retain full responsibility for appropriate management of any situation. Sonicaid Fetalcare 3 is designed as an information system intended to present information to assist clinicians in delivering the highest possible standard of care, not to replace established clinical practice. All users are responsible for ensuring the accuracy of entered data, and for confirming that it has been correctly logged.

**Date / Time**

All actions, traces, data entry, etc. are time stamped using the server system clock. If the system clock is incorrectly set, the logged times will reflect this error. The user is responsible for checking that the date & time are correct – this is shown in the bottom right-hand corner of the screen at all times. If the time is incorrect, advise the system administrator immediately.

**Patient details**

Any patient names, demographic or other data appearing in this document are for demonstration purposes only and are purely fictitious. Any similarity this data may have to any real person is purely coincidental.

**System operation**

The system is designed for continuous operation and should therefore never be switched off during normal use.

**CTG analysis (option)**

This document only describes the operation of the analysis function. It does not include the clinical application which is beyond the scope of this document. It is essential that all users of the analysis option are fully trained on its use and application. This function is an aid to CTG trace interpretation to inform the clinician. It does not diagnose fetal condition and does not replace the need for expert trace interpretation & effective management of the pregnancy. It is only approved for use from 26 weeks to term, prior to the onset of labour, and is not for use during labour. Visit www.huntleigh-diagnostics.com for further information on the CTG analysis.

**Potential loss of clinical data**


Installing FC3 on laptops or other mobile devices is not recommended. If used, power saving actions relating to lid closure, & inactivity timeout functions (standby / sleep / hibernation) must be disabled & mains power must be used at all times. **Failure to disable such features, or battery operation, may result in irretrievable loss of clinical data.** The use of laptops or other mobile devices as extra client/viewing terminals is acceptable. However, it is recommended that the above features are disabled on clients as well, to ensure continuity of operation. **Desktop PCs must also have any sleep / hibernation / standby functions disabled to avoid data loss.**

The use of screen savers is acceptable. However, users must be aware that when these are active, update of live data is suspended. Screen savers which superimpose on, but do not obscure, the trace view (eg. MS Windows “Bubbles”) may cause confusion as the FC3 view can still be seen but will not be updated. However, no data is lost in screen saver mode & the view will be updated as soon as the screen saver is deactivated.

2 Getting started

2.1 Screen layout

The screen is divided into two areas, as shown below:



- The left hand view displays bed group(s), listing the beds contained within the system.
- The right hand view will contain one of the following at any one time
 - Home page
 - Single trace view
 - Multi trace view
 - Patient record view
 - Analysis results view (where option installed)
 - Analysis trend view (where option installed)
 - Chalkboard (where option installed)
 - Remote CTG (where option installed)
 - Security and audit (where option installed)
 - GDT (where option installed)
 - Partogram (where option installed)

3 System operation

Starting the Sonicaid Fetalcare 3 application will show the ‘Home page’ as displayed below.


<h3>3.1 Home page</h3>	
	<p>Note: this screen may vary in different markets</p>
<h4>3.1.1 To return to the home page</h4>	
	<p>Click to return to the home page.</p>
<h4>3.1.2 Help</h4>	
	<p>Click for ‘Instructions for use’</p>
<h3>Bed frame view</h3>	
<h3>3.2 Bed & bed group functions</h3>	
	<p>The left-hand column shows all the beds in the</p>

	<p>system, organised into bed groups.</p> <p>Click on a bed to admit a patient or view a CTG trace.</p>
	<p>The bed group shows the group name and number of occupied beds / total beds in the group.</p> <p>Click to maximise or minimise the beds in this group. This turns red when any alert is active in this bed group.</p>

3.3 Adding a new patient:



Select a bed.



Enter 'Last name', 'First name' and 'Hospital ID'.



As the patient details do not match an existing patient within the system, a 'Create new' patient button appears. Pressing this button adds the new patient details to the system and displays the chalkboard entry view.

Enter Chalkboard Data

Consultant

Comments

Midwife

Details can be added to the chalkboard at this point and/or at a later date.

Antenatal 2/4

- Bay 1
- Bay 2
- Bay 3 Jones, Mary ♥ 135
- Bay 4 Williams, Brenda ♥ 130 ♥ 139

Labour 1/4

- Room 1
- Room 2
- Room 3
- Room 4 Foster, Mary ♥ 127 ♥ 151 ♥ 121

Theatre 4/4

- Theatre 1 Pearson, Mary ♥ 155
- Theatre 2 Davis, Mavis ♥ 155
- Theatre 3 Wright, Mila ♥ 156

Current Pregnancy

Pregnancy Menu | Pregnancy Notes | CTG 29/01/2016 10:54

♥ 155 UA 21

200 180 160 140 120 100 80 60 US

10:54 (29/01/2016) 1cm/min 11:04 (29/01/2016) 100%

Toco (%)

The patient is booked into the selected bed.

3.4 Admit an existing patient



Select a bed.

Select Patient

Enter patient details

Last Name

First Name

Hospital ID

NHS ID


Matching Patient List

Last Name	First Name	Hospital ID	NHS ID	Date of birth	Admitted
Smith	Lisa	93847545			
Smith	Mary	82743654			
Smith	Sally	47839781			



Enter patient details, matching patients will be displayed as you type.

Then click on the correct patient name in the 'Matching Patient List' and click on the '✓' button, or double click on the patient name.


IMPORTANT:
Always check that the patient I.D. is correct – this is the only unique identifier, as there may be more than one patient in the database with the same name.



The bed symbol changes from grey to blue, indicating that the patient is now admitted to the bed.





Bed frame showing bed name, patient name and the heart rate data from the CTG. The bed frame border colour changes to green to show that this is the currently selected bed.

 <p>Bay 4 Williams, Brenda 1 142 2 154</p>	<p>Bed frame showing twin and triplets heart rates.</p>
 <p>Room 4 Foster, Mary 1 157 2 147 3 109</p>	

3.5 Multi-bed view

	<p>Pressing this button displays the multi-bed view settings screen</p>
---	---

 <p>Multi Bed View Settings</p> <p>Number of columns</p> <p><input type="radio"/> 1 <input checked="" type="radio"/> 2</p> <p>Select bed groups to display</p> <ul style="list-style-type: none"> <input checked="" type="radio"/> All traces <input type="radio"/> Custom <input type="radio"/> Antenatal <input type="radio"/> Labour <input type="radio"/> Theatre <p style="text-align: right;"> <input type="button" value="X"/> <input type="button" value="✓"/> </p>	<p>Select the number of columns view required (1 or 2).</p> <p>Select the bed groups to display. This section allows the user to select specific beds to display in the multi bed view window. There are three options available.</p> <p>1 - The user can select all traces to display all beds with live CTG traces associated.</p> <p>2 – The user can select a specific bed group. All beds with live CTG traces associated with them within the selected group will be displayed.</p> <p>3 – The user can select ‘custom’. A mixture of beds</p>
---	--

	<p>with live CTG traces associated with them can be selected from across all the bed groups.</p> <p>The multi-bed view displays all beds with active traces.</p> <p>The icons in the header bar for each trace operate in the same way as in single trace view.</p> <p>To select the single trace view for one of the traces, double click on the trace.</p>
---	--

3.6 Reviewing stored patient data

	<p>Review stored patient data for a patient not currently booked into any bed.</p>
---	--

3.7 Starting a CTG trace



Check that the fetal monitor is connected to the wall connection point.

Start monitoring the patient. This will automatically start the trace recording session.



If a trace is completed without admitting a patient, the trace will be unallocated. Contact your system administrator.

3.8 Stopping a CTG trace



Switch the fetal monitor off. After a short delay the CTG will close.


 If the patient has left the bed, then you must remember to discharge the patient on the system.



If a new trace is started on a new patient,

without first discharging the previous patient, the trace will be assigned to the previous patient. Contact your system administrator.

3.9 Single trace view



This is the default single patient view.



Tab shows a patient name.



Shows bed name and FHR data.

3.9.1 Trace details

The chart speed can be set to 1, 2 or 3cm/min.

FHR vertical scaling can be set to 20bpm/cm, over the range 50-210bpm, or 30bpm, range 30-240bpm.

TOCO scale covers the range 0 to 100% for external TOCO and 0 to 100mmHg for IUP monitoring.

When available, maternal heart rate (MHR) will be plotted in green on the FHR scale. All other maternal parameters will be displayed on the contractions scale.







Chart date & time are displayed at fixed intervals between the FHR & TOCO scales. The interval is 10 minutes when the chart speed is set to 1cm/min.




The FHR mode is annotated towards the bottom of the FHR chart as 'US' or 'FECG'.





3.9.2 Trace scrolling



3.9.3 Alerts

Note that these alerts operate independently of the local alert system of each connected fetal monitor and the settings may be different.



These alerts are not clinical alerts, they are user alerts.

Although it has become standard practice in the industry to refer to these as, for example ‘Tachycardia alarm’, implying by its name some clinical significance, this is not the case. They do not interpret the FHR data in any way, they simply draw the user’s attention to the fact that the FHR has been outside a user defined range for a user defined period of time.

Additionally, the Low & High FHR alerts make limited allowance for loss of signal & transient returns to rates within the user set limits. However, the presence of signal loss or transient returns to rates within the user set limits may result in these alerts not triggering.


In the event of an alert, it remains the responsibility of the user to determine what has caused the alert, whether there is any clinical risk, and to ensure appropriate management.

Alerts can be disabled by the user, or the audio can be turned down or off.

The alerts in this product are provided as an adjunct to good clinical practice in checking the patient & reviewing the trace on a regular basis & must not be relied on to detect abnormalities in the trace.



Click to check or change alert settings.



The system includes alerts for the following:



High FHR alert - alerts to FHR above the set threshold for a set time period.



Low FHR alert - alerts to FHR below the set threshold for a set time period.



Loss of signal alert - alerts to loss of FHR signal over a set percentage of a set time period.

The system also includes a cross-channel alert. This alerts users to conditions where the heart rate is the same (or similar) on any 2 heart rate channels (FHR1, FHR2, MHR).

Note that the cross-channel alert cannot be adjusted and is permanently enabled.



Alert threshold settings are shown on the chart by a change in background colour.

Alerts are indicated on the chart:




Indicates an alert event.




Indicates alert event acknowledged

Hover over the

	<p>icon to show further details of the alert.</p> <p>If an alert is active on a bed, the heart symbol and the FHR data will flash</p>
	<p>Red bell indicates active alert, click on it to acknowledge the alert.</p>
<p>3.9.4 Adding a trace annotation</p>	
	<p>Click to add a CTG annotation.</p>
	<p>There are 3 data fields:</p> <ol style="list-style-type: none"> Event: Click on the drop-down box for a pre-defined list of annotations. <p>Alternatively, simply type free text into this field.</p>


	<ol style="list-style-type: none"> 2. Details: Type details into this field. 3. Author: Enter your name.
	<p>An annotation is shown on the CTG. If this is added more than 1 hour after the event, or on a stored trace, the mark will be grey instead of yellow.</p> <p>Hover over the annotation to view the details.</p>
<h3>3.9.5 Printing a trace</h3>	
	<p>Click to print the CTG.</p> <p>The last page of a print-out shows a list of all annotation details.</p>

3.9.6 Connecting SRX series Doppler - optional



If non-rechargeable batteries are being used on the device, check the Doppler has had the battery type configured appropriately so as not to charge a non-rechargeable battery

Connect the SRX Doppler using a micro USB cable to the USB connection on the workstation. The first time this is connected the device driver should be loaded.



Switch on the Doppler, the 'SR Doppler' icon should appear showing the serial number at the top of the left hand side.

	<p>Click on the 'SR Doppler' button to access stored and real time trace</p>
	<p>Click 'Doppler Record' to see real time data or the icon to the right to allocate the recording to a patient entering patient data</p>
	<p>Click 'Associate this session with a bed' if a patient you are recording is already in a bed but does not have a CTG connected. Click 'Reset this session' to start new recording</p>

<p>10:35 Doppler Record  Stored. USB. Started 28/11/2018 10:35 Not assigned to a patient.</p> <p>10:38 Upload all stored records by Administrator</p> <p>10:38 Doppler Record  Stored. Uploaded. Started 20/10/2016 07:09 Not assigned to a patient.</p> <p>10:38 Doppler Record  Live. USB. Started 28/11/2018 10:38 Not assigned to a patient.</p>	<p>Select 'upload all stored records', click 'Doppler Record' to open recording and the icon to the right to allocate the recording to a patient</p>
---	--

<p>Select recent session to view</p> <ul style="list-style-type: none"> Current session: 28/11/2018 10:35 Current session: 28/11/2018 10:35 Previous session: 28/11/2018 10:28 Previous session: 27/11/2018 16:41 Previous session: 27/11/2018 16:39 Previous session: 27/11/2018 16:35 	<p>From the drop down box select the connection session</p>
---	---

<p style="text-align: center;"><u>Details</u> </p> <p>Monitor Serial Number: C200573917</p> <p>ID Code: SRX</p> <p>First connected: 11/27/2018 4:35:37 PM</p> <p>Display Name: C200573917</p> <p>Description: <input type="text"/></p>	<p>Details of the Doppler should be shown, clicking the 'Details' icon will allow you to change the name displayed and enter a description for the device. This is stored against the device ID on the system</p>
--	---



3.9.7 CTG analysis (optional feature)





CTG Analysis Intended Use

The intended use of the CTG analysis option is for the computerised analysis of antepartum cardiotocograms in pregnancies from 26 weeks gestation onwards (32 weeks in the USA). It can be used on women who are experiencing Braxton-Hicks contractions but is not intended for use in established labour as the fetus is then exposed to additional factors such as labour contractions, pharmacological agents, and epidural anaesthesia.

The analysis provided is intended as an adjunct to - and not a replacement for - the clinician's visual assessment of a cardiotocogram. As such, this CTG analysis is an aid to clinical management but not a diagnosis, which remains the responsibility of an appropriately qualified clinician. Indeed, both the clinician's visual assessment of a cardiotocogram and the analysis provided by this software should be considered within the context of a full clinical assessment before decisions are made regarding management. Such an assessment may include further tests such as umbilical blood flow velocity waveforms or biophysical profiling.


	<p>Click to start a CTG analysis.</p>
	<p>Enter the gestational age (GA).</p> <p>Select 'One Mark' to start a CTG analysis.</p> <p>Select 'Two Marks' to start a CTG analysis for a defined section of the CTG Minimum period 10 minutes, maximum 60.</p> <p>Click '✓' to confirm that the patient is NOT in established labour.</p>
	<p>Should you choose to proceed with a CTG analysis while the patient is in labour, you MUST NOT rely on the analysis result.</p>

	<p>For 'One Mark' analysis, place the vertical marker at the desired start point and click.</p> <p>Note that there may be up to 10 minutes delay before any analysis results become available.</p>
---	--

	<p>For 'Two Mark' analysis, place the first vertical marker at the desired start point and click. Move to the desired end point and click to place the second marker, which will move in 2 minute steps.</p>
---	--

	<p>Indicates the start point of the analysis.</p>
---	---

	<p>Indicates the end point of the analysis.</p>
---	---



FHR1
High Episodes
Low Episodes
Signal loss >25%

Latest 01/02/2016 11:01:48

Started 01/02/2016 10:49:48

ADVIC

Next result (01:19)

Results	12	10
Signal loss (%)	0	0
Contractions	0	0
Fetal movements (per hour)	0	0
Moves per min high	0.0	0.0
Moves per min low	0.0	0.0
Basal Heart Rate (bpm)	154	136
Accelerations > 10 (bpm)	0	2
All decelerations	2	1
Significant decelerations >20 (lost beats)	1	0
Area of largest deceleration (lost beats)	28	10
Minor decelerations <=20 (lost beats)	1	1
High episodes (mins)	0	0
Low episodes (mins)	0	0
STV overall (ms)	4.6	5.1

If analysis results are not ready the screen will show “Awaiting results”. The first result will be displayed after 10 minutes of good quality trace. This is updated every 2 minutes up to a maximum of 60 minutes.

The screen shows a compressed CTG with analysis results below.

Dawes Redman Criteria MET by FHR1 at 14 minutes

Dawes Redman Criteria NOT YET MET by FHR1 at 12 minutes

Dawes Redman Criteria NOT MET by FHR1 at 24 minutes

There are 3 possible outcomes as per these examples. Refer to the clinical training material for further information on understanding these outcomes.

Criteria not met because:

- No moves and less than 3 accelerations
- Baseline fitting is uncertain

Where criteria is not met, or not yet met, moving the mouse pointer over the coloured result bar (above) will show the reasons for the criteria not being met.



For twins analyses click the tabs to switch between each fetus' analysis.


The tabs are colour coded in the same way as the results, indicating whether the latest result meets the criteria or not.




Click to print the analysis.



Click to stop a live analysis.



Click to select the trend view.




The trend view presents graphic plots of 3 key analysis parameters against time:

- STV
- Fetal Movements
- Basal heart rate

A minimum of 3 sets of traces and analyses are required to activate this function.

Note: Trends are not available for twins.



Click to select time period for the trend view.



Click to print the trend view.

3.9.8 STAN



STAN monitoring

The STAN functionality is an aid to clinical management but not a diagnosis, which remains the responsibility of an appropriately qualified clinician. Sonicaid Fetalcare 3 is merely presenting data from the STAN monitor, and therefore, any clinical significance of the data is down to the user, who should refer to the STAN documentation and Neoventa's guidelines.





It is essential that users of this system are fully trained in using ST analysis and in understanding and interpreting such data. Huntleigh accept no responsibility for this & users are referred to Neoventa for all support in relation to this analysis.






This option allows Sonicaid Fetalcare 3 to receive ST analysis data from Neovanta STAN® fetal monitors when connected.

This data is displayed with the CTG trace & is archived as part of the CTG trace record.


NOTE: The scaling of the STAN area of the trace screen is different from that of a STAN monitor. Users who are familiar with STAN should be aware that data will be scaled differently on Sonicaid Fetalcare 3 and could give the impression that the ST rise is lower than it actually is.



The t/qrs ratio data associated with the ST analysis is displayed on a separate chart below the contractions chart in the CTG trace view.


	<p>ST events are displayed at the top of the FHR scale.</p> <p>Hover over the ST event mark to display additional information relating to the ST event.</p>
<h3>3.10 Switching views</h3>	
	<p>Click to select CTG view.</p>
	<p>Click to select a stored CTG.</p>
	<p>Click to view patient details.</p>
	<p>Click to view analysis results</p>
	<p>Click to view analysis trends (only active when 3 or more trace analyses completed)</p>

3.11 Patient details view



This view shows all activity associated with the patient.


3.11.1 Editing patient details



Press the edit details button.


<p> Edit Patient Details</p> <p>Last Name <input type="text" value="Brook"/></p> <p>First Name <input type="text" value="Emma"/></p> <p>Maiden Name <input type="text"/></p> <p>Hospital ID <input type="text" value="93248752"/></p> <p>NHS ID <input type="text"/></p> <p>Date of birth <input type="text" value="DD"/> <input type="text" value="MM"/> <input type="text" value="YYYY"/></p> <p>Doctor <input type="text"/></p> <p>Address 1 <input type="text"/></p> <p>Address 2 <input type="text"/></p> <p>Address 3 <input type="text"/></p> <p>Address 4 <input type="text"/></p> <p>Tel 1 <input type="text"/></p> <p>Tel 2 <input type="text"/></p> <p></p>	<p>Edit details as required.</p>
---	----------------------------------

3.11.2 Retrieving stored CTGs





Click on the trace to be retrieved.

3.11.3 Retrieving stored analysis results



Click to expand the CTG details and click on the analysis to be reviewed.

3.11.4 Adding patient notes

	<p>To add patient notes</p>
	<p>Type notes as required</p>

3.12 Chalkboard



Click to open the chalkboard.

Chalkboard	Consultant	Comments	Midwife
Antenatal 2/4			
Bay 1 Evans, Phoebe			
Bay 2			
Bay 3 Jones, Mandy			
Bay 4			
Labour 0/4			
Room 1			
Room 2			
Room 3			
Room 4			
Theatre 4/4			
Theatre 1 Clark, Maria			
Theatre 2 Jones, Mat			
Theatre 3 Davis, Sally			

The chalkboard provides an overview of the current status of all beds. Each bed has a chalkboard entry associated against it. To add or edit chalkboard data, click on the desired chalkboard field.

Edit chalkboard data

Consultant

Comments

Midwife

Enter or edit data in the chalkboard data view.

If a bed has no patient associated to it, chalkboard data cannot be added to that entry.

3.13 Discharging a patient

Click to discharge a patient from a bed.

3.14 Settings, Audit and Admin

The settings available to change for a user will depend on the user access level.



Select the 'settings, audit and admin' button.

Four buttons appear:

Admin access.
User settings.
System audit.
GDT settings.



System Administration

Current Users | Users | Machines | System Access Levels | System Settings | Bed Group Profiles | Log Messages | Features

User Logon Id	User	Licence Id	Logged On At	Machine Id
9	Administrator	0	29/02/2016 08:55	2

View User | View Machine

Admin access.

Current user:
A list of current users logged on to the system.

System Administration

Current Users | Users | Machines | System Access Levels | System Settings | Bed Group Profiles | Log Messages | Features

User Id	User Name	Display Name	System Access Level Id	Confidentiality Level	Blocked	Account Inactive
1	Centrale	Centrale	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
2	AutoDischarge	AutoDischarge	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
3	NoSecurity	NoSecurity	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
4	Monitor	Monitor	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
5	STAN	STAN	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
6	Reserved6	Reserved6	Administrator	High	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
7	Reserved7	Reserved7	Administrator	High	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
8	Reserved8	Reserved8	Administrator	High	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
9	Reserved9	Reserved9	Administrator	High	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
10	Reserved10	Reserved10	Administrator	High	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
11	Administrator	Administrator	Administrator	High	<input checked="" type="checkbox"/>	<input type="checkbox"/>
12	matt.clarke	Administrator	Administrator	Low	<input type="checkbox"/>	<input type="checkbox"/>

Create New | Edit | View System Access

Users: A list of current users and their user access levels.

System Administration

Current Users | Users | Machines | System Access Levels | System Settings | Bed Group Profiles | Log Messages | Features

Machine Id	Machine Name	System Access Level Id	Confidentiality Level	Static Licence	Blocked	Deleted
1	CentraleServer	Administrator	High	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
2	I10-1764	Administrator	Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Create New | Edit | View System Access

Machines: A list of machines which have access to the system.

System access levels: Create and edit system access levels.

The screenshots show the following sections of the System Administration interface:

- System Access Levels:** A table with columns 'Name' and 'Description'.

Name	Description
Administrator	Full user access, full admin access
Midwife	Full user access, no admin access
- CTG Settings:** A form with the following settings:
 - Chart Speed: 1cm/min
 - FHR Vertical Scale: 20 bpm/cm (50-210 bpm)
 - IUP Scale: mmHg
 - Minimum Graph Width (cms): 30
- Bed Group Profiles:** A table with columns 'Name' and 'Description'.

Name	Description
All Bed Groups	All Bed Groups
- Log Messages:** A table titled 'All Log Sessions' with columns 'Id', 'Source App', 'Session Start', and 'Is Live'.

Id	Source App	Session Start	Is Live
6	CentralClient_110-1764	23/02/2016 15:05	●
7	CentralClient_110-1764	24/02/2016 10:45	●
8	CIS	25/02/2016 08:21	●

System settings:
System wide CTG settings and general settings.

Bed group profile: Create new and edit existing bed group profiles.

Log messages: Displays any errors that have occurred.

Features: Displays all features available and which ones have been licensed.

	
	<p>User settings.</p> <p>Allows the logged on user to set specific CTG settings and which bed group profile they wish to use (if any).</p>
	<p>System audit.</p> <p>A comprehensive audit of all user transactions.</p>
	<p>GDT settings.</p> <p>Enable GDT. Setup incoming and outgoing GDT file locations. Select which file format extension is to be used.</p>

4 Trouble-shooting

Due to the nature of the system, it is not possible to cover all possible areas of trouble-shooting in this manual. This section includes first-line trouble-shooting relating primarily to difficulties in using the system. For more in-depth system support refer to your system administrator.

Problem	Possible solutions / explanations
Screen is blank	<ul style="list-style-type: none"> -Check screen connected to mains power. -Check power is switched on at mains socket -Check screen on/off switch is On. On most screens there is a small green or amber light – confirm that this is lit. -Blank screensaver may have been activated – move/click mouse or enter keystroke on keyboard to re-activate screen -Computer may have been switched off.
Not receiving data from CTG	<ul style="list-style-type: none"> -Check CTG is switched on and working (note – CTG printer does not need to be running) -Check CTG is connected to wall socket in room – check cable connectors are secure -Book patient into bed – traces cannot be viewed until this is done
Unable to annotate trace	<ul style="list-style-type: none"> -Traces can only be annotated within the timeframe of the trace, not on the blank grid area to the right of the trace.

System not responding	<p>Check connections between the computer, keyboard, mouse and the network socket.</p> <p>Refer to your administrator – the server may need to be shut-down and re-started – do NOT attempt to do this unless authorised and trained.</p>
All access terminals shut-down and not responding	<p>Power cut – the main server will be supported for a short time by the UPS (depending on model) – typically about 10-15 minutes. After this, if power is not restored, the server will be shut-down. When power is subsequently restored, the whole system will need to be re-started – contact your administrator or IT department</p> <p>If power is restored before the server shuts down, simply re-start all access terminals.</p>
Mouse / Keyboard not responding	<p>Check cable and connection</p>
Print-outs not printing	<p>Check printer is switched on, is ‘On-line’ and has sufficient paper loaded.</p> <p>For ink-jet printers, ink cartridge may need replacing – refer to printer manual for details.</p> <p>For laser printers, toner may need replacing – refer to printer manual for details.</p>
System error message appears on screen	<p>During system maintenance, while shutting down or re-starting the system, error messages may appear – these will normally clear themselves after a short delay. If the message does not clear after ~1 minute, contact your system administrator or first line support team.</p>

5 System maintenance

5.1 General maintenance

The only scheduled maintenance required is to back up the database. This is subject to local policy and is normally managed by your IT department. Contact your system administrator or IT department for details.

IMPORTANT: In the event of hardware failure, software bugs or other system related problems, disk storage overflow, etc., data may be lost at any time. Note that data may also be lost in the event of network issues and other infrastructure issues which are managed by and remain the responsibility of the customer. The customer is responsible for ensuring that regular back-ups of the database are kept in accordance with established industry practice, local protocols and guidelines. Failure to do so may result in total loss of all patient information, CTG trace records, etc.

Huntleigh Healthcare cannot, under any circumstances, accept any responsibility for loss of, or corruption of, any stored data relating to Sonicaid Fetalcare 3. Such data is the property of the customer who is solely responsible for protecting it.

No other maintenance is required, other than cleaning and checking for damage. This can be included in your routine local equipment maintenance programmes.

If you have maintenance contract cover, contact your system administrator or supplier for details. Note that this does not include data back-up which at all times remains the responsibility of the customer.

6 Service support

First-line support is provided by your local super-user or local IT support. For service support & details of maintenance contracts, contact your supplier.

6.1 Licence renewal

The use of this software is controlled by a licence key. To extend or renew your licence, or to add software options, refer to the installation instructions supplied with the product. Contact your supplier for further information on options, upgrades & support.



Medical Devices Directive 93/42/EEC

MANUFACTURED & DISTRIBUTED IN THE UK BY:

*Huntleigh Healthcare Ltd
Diagnostic Products Division
Cardiff CF24 5HN UK
Tel: +44 (0)29 20485885
Fax: +44 (0)29 20492520*

*Email: sales@huntleigh-diagnostics.co.uk
www.huntleigh-diagnostics.com*

As part of the ongoing development programme, the company reserves the right to modify specifications and materials of the Sonicaid Fetalcare 3 without notice.

Huntleigh Healthcare Ltd – A member of the Getinge group. Sonicaid, Huntleigh and ‘H’ logo are registered trademarks of Huntleigh Technology Ltd.

©Huntleigh Healthcare Ltd. 2016.